


# Leads for adhesive electrodes with snap fastener



*Example of use*



- Simple reconnection of the lines at the distributor head of the support arm (by replacing the suction lines)
- Suitable for use with all our suction systems
- For use e.g. in infectious patients or children
- No reconnection at the ECG device necessary due to temporary use of a patient cable
- No need to sacrifice the convenience of a support arm

Length of leads C1-C6 = each 1,2m

Length of leads R, L, F, N = each 1,4m

Available both individually and as a set (10 lines)

Only suitable for use with adhesive electrodes with snap-fastener!

strässle&co



Medizintechnik GmbH

**Electrode cables for**

**adhesive electrodes with snap-fastener**

Dear customers,

are you aware that our suction systems can be supplemented by interchangeable **electrode cables for use with adhesive electrodes?**

With the aid of these electrode cables, it is possible to temporarily switch to adhesive electrodes for use on high-risk patients or children by simply reconnecting the cables to the distributor head of our support arms, without having to forego the practical conveniences of a support arm. This also saves the user having to switch to a patient cable in the meantime, and there is no need to reconnect the ECG device. This means that users can continue to use their familiar equipment.

These electrode cables can be used on all our existing suction systems, regardless of the type and age of the equipment. The electrode cables are only suitable for use with commercially available snap-fastener electrodes.

**Please find further informations on back side.**